



## CHEMSTRESS®: Biomanufacturing Process Control

Chinese hamster ovary (CHO) cells, the powerhouse of the biopharmaceutical industry, are inherently variable. This functional heterogeneity is what makes them such a valuable tool for producing biotherapeutics, while also accounting for the considerable cost and time associated with cell line development (CLD). As is the case with other transformed mammalian cell lines maintained *in vitro*, inherent genomic instability can modulate chromosomal arrangement and complement gene copy number and transcriptional activity. In industry, exploitation of clonal genetic variation is essential to derive cell lines capable of supporting biomanufacturing, as this underpins a variety of prerequisites, most particularly an ability to survive in synthetic environments, achieve high rates of cell proliferation and efficiently manufacture recombinant proteins<sup>1</sup>. However, genetic variation can also give rise to unpredictable behaviour *in vitro* such as loss of productivity<sup>2</sup>, variations in product modifications, e.g. *N*-glycosylation<sup>3</sup>, and proliferation rate<sup>4</sup>. Dealing with the consequences of clone-specific genetic and functional variation<sup>5</sup>, and its acquisition during CLD, can be extremely time-consuming and labour-intensive.

CLD requires a significant amount of investment and the instability of cells is the central driver of this lengthy process. It is essential to have the appropriate instrumentation, expertise and time, to assess the performance attributes and stability of hundreds of clonally-derived cell populations in parallel. In addition, it is necessary to culture these cells in an optimum cell environment, to maintain a stable state. Development of cell culture media is therefore another major source of continuing investment. This involves the screening of chemical additives, the functional validation of these enhancer molecules, and the confirmation of the production of a stable, robust, high-quality media product. Moreover, quality control (QC) of the complex medium product itself (post-production, transport and storage) is generally restricted to simple monitoring of cell growth and productivity. There is an urgent requirement for new QC technologies able to measure lot-to-lot variation and functional integrity of this vital part of the bioproduction process.

As the biopharmaceutical market expands to allow for more competitive conditions in the form of biosimilars of costly blockbuster drugs, new technologies which can decrease the time and costs associated with the development of these biotherapeutics will be essential. These technologies must be prepared to understand the molecular mechanisms of the cells at the centre of this industry. To illustrate the financial consequences of defects in any one of the steps in the bioproduction process, contamination alone can

result in an estimated \$1bn loss in revenue for blockbuster biotherapeutics, along with a projected \$1m loss in production costs and up to \$20,000 in quality assurance investigations<sup>6</sup>. Companies are therefore ultimately driven towards streamlining CLD and the bioprocessing of biotherapeutics, through the use of informative assays which can generate functional data and determine reliable QC of reagents.

### ChemStress®: Generating a Functionally Relevant Biological Signature

ChemStress® is an information-rich, analytical assay supplying data on the functional quality of clones, facilitating the selection of stable producer cells and high-quality media environments through the generation of cellular responses to chemical stressors. The technology draws on the unique biological signatures (e.g. growth and titer for CHO cells producing monoclonal antibodies) a given cell generates when it is grown in the presence of small molecule stressors, which have been selected to simulate toxic environments or stimulate specific stress response pathways. This data gives biological information on how specific clones respond to stressors, and whether these clones have the ability to remain stable and productive in these conditions. These cellular responses are unique identifiers, or 'fingerprints'. While classic fingerprints are static measurements of identity, these 'ChemStress® fingerprints' are markers which, much like cells, can change overtime as a result of e.g. epigenetic silencing events, or loss of recombinant protein expression, giving early indications of loss of specific functionality. The diverse range of responses is a unique way of: identifying cell populations, indicating divergence from its original state and indicating the growth quality of cell culture conditions. ChemStress® employs the use of a novel 96-well microtiter plate, wherein each well is coated with a single specific chemical (Fig. 1), selected based on its published ability to simulate specific conditions, or effect specific cellular pathways. The workflow is illustrated in Fig. 2, where product titer is analysed using our Valita®TITER assay<sup>7</sup>.





Figure 1. ChemStress® technology. Each well of the 96-well plate is coated with a single specific small molecule stressor, at select concentrations. Each chemical is selected to stimulate a particular cellular pathway, e.g. apoptosis and oxidative stress.



**Figure 2. ChemStress® assay.** Cells are seeded into the plate at a particular density, depending on the cell system being examined. Following three days of cell culture, the cells' growth and titer response are determined from the contents of each individual well and compared to control wells without stressors. In the current example, a CHO-cell system producing monoclonal antibodies (mAbs), the manufacturing metrics are growth and product titer. Following analysis of the data using proprietary software (ValitaApp®), a ChemStress® fingerprint is generated.

The assay was developed in collaboration with researchers from the University of Sheffield and chemicals were selected based on their published effect on cellular systems (examples included in Fig. 1). Chemicals were tested against a mAb-producing CHO cell line and screened using half maximal inhibitory concentration ( $IC_{50}$ ) assays to determine cell growth inhibition. Several chemicals are present at multiple concentrations to allow for capture of variable responses when exposed to different concentrations of key chemicals.

The wide range of selected chemicals are selected to target a number of pathways stimulated under stressful conditions, including nutrient starvation, cell death and oxidative stress. Rapamycin, for example, is an inhibitor of mammalian target of rapamycin, which is suppressed during nutrient starvation<sup>34</sup>. Rapamycin induces autophagy, a pathway which causes controlled breakdown of cellular compartments, or removal of damaged organelles, generating the components required by the cell for survival during nutrient starvation<sup>35</sup>. Brefeldin A is often used in flow cytometry to inhibit protein secretion and accumulate cytokines for staining<sup>14,15</sup>. It is also used to induce endoplasmic reticulum (ER) stress and eventual apoptosis<sup>16</sup>, in addition to causing disassembly of the Golgi and fusion with the ER<sup>17</sup>. There are a number of other chemicals on the ChemStress® plate which induce apoptosis, including cadmium acetate and menadione sodium bisulfate, and chemicals which induce oxidative stress, e.g. cobalt chloride and 3-amino-1,2,4-triazole.

## Applications


ChemStress® targets core issues in biomanufacturing, including:

- Identify parental cell line of clonal isolates (Valita®ID),
- Stability of producer cells (Valita®STABILITY),
- Quality control of cell culture environments (Valita®QC).

## Valita®ID


Cross-contamination of cell lines is a major issue, not just for the biopharmaceutical industry, but research in general. HeLa cell contamination of multiple cell lines is one such well-documented case. For the biopharmaceutical industry, regulators insist that the identity of clonal isolates can be traced. There are techniques available to allow companies to meet these requirements, which include isoenzyme analysis, karyotyping, short tandem repeat (STR) analysis and DNA sequencing and barcoding. ChemStress® fingerprinting exploits the functional variability of cell populations to

generate biologically relevant data on the specific clonal identity of cells, giving users a quick and easy method to keep track of cells (illustrated in Fig. 3).



**Figure 3. Valita®ID workflow.** Prior to banking of cell lines, ChemStress® fingerprints are generated and stored in a database. When cells are used during the CLD process, it is possible to generate a fingerprint of a query cell to determine the identity of the cell line from which the clonal isolate derives.

To test this, ten different mAb-producing CHO cell lines were sourced, derived from a variety of different CHO cell families, including CHO-K1, CHO-S and CHO-DG44. Each cell line was coded anonymously as CHO-A to CHO-J, to prevent any bias. ChemStress® profiling was able to cluster the different groups of coded cell types together as indicated in Fig. 4. The analysis used 'linear discriminant analysis (LDA)' and 'Euclidean distance minimisation' in addition to 'Ward's Hierarchical clustering' and found that fingerprints were unique and could be used effectively to partition cell types into distinct parental 'groups'. Also, using LDA, given a query fingerprint, the probability of it being from a specific group can be approximated. This application can be utilised to generate a library of clonal fingerprints to allow for the identification of clones following storage and during CLD.



**Figure 4. Cluster Dendrogram** shows the correct clustering of the cell lines based on the ChemStress® technology. It is possible to determine the parent cell line of specific clonal isolates. This allows users to keep detailed biological records on clone identity and to trace clones back to the specific isolate from which they were generated.

Given the need to meet the regulatory requirements on clonal identity, it is necessary for the generation of easy to implement, simple, rapid and cost-effective analytical assays which can determine the origin of clones. The ChemStress® platform offers a simple assay which generates biologically relevant data, that predicts the parental identity of clonal isolates.

## Valita®QC

There is a substantial lack of informative, cost-effective assays to allow for screening of cell culture medium prior to usage and following long-term on-site storage. Lot-to-lot variations in media have been well-documented, in addition to the consequences of incorrect storage including temperature fluctuations and UV-exposure<sup>36,37</sup>. Valita®QC can be used to



screen the growth and productivity of a reference cell line when grown in the presence of query lots of media, generating informative data which allows for the acceptance or rejection of media lots based on functionally relevant responses (concept illustrated in Fig. 5).




Figure 5. Valita® QC workflow. A cell line is grown in different lots of a specific media. The resulting fingerprints are analysed to determine acceptable responses (acceptable growth and productivity) to allow for the acceptance of the lots; those with poor fingerprints are rejected.




Figure 6. Growth responses following culturing of cells in different media lots. Three growth response graphs are indicated from three specific chemicals on the ChemStress® plate. The same cells are shown to have unique profiles, despite being cultured in different lots of the same composite media, clearly indicating quality difference between the lots.

The growth and titer responses of cells cultured in different lots results in a distinct fingerprint, as illustrated from select chemicals in Fig. 6. The overall fingerprint generated from the multiple chemicals screened can be used to determine if lots should be rejected based on poor growth and productivity of cells. In addition to determining lot-to-lot variation, ChemStress® aims to detect if any damage has occurred to affect the quality of the media prior to use, e.g. following incorrect storage, including exposure to UV or temperature fluctuations, or during long-term storage (ageing). It is important to have a functionally relevant, reliable assessment of the quality of media before the product is released or used. As a critical component of the end stages of a lengthy process, it is necessary to have stringent QC technology in place which can allow for the accurate assessment of this essential reagent.

## Valita® STABILITY

Predictable bioprocesses require stable producer cells that are not prone to sudden changes leading to diminished productivity or quality. To demonstrate stability, clones are sub-cultured and monitored for titer, growth and various CQAs over a period of months. There is a very high workload associated with these trials. To enhance trial efficiency, ChemStress® is being used to identify maximally stable clones pre-trial and to fail unstable clones quickly within a trial.

As a broad measure of stability, companies typically monitor titer over a period roughly equivalent to a production run, with any clone losing more than 30% of its initial titer over this time being considered unstable. Maintaining titer can be a misleading stability metric if a clone's cell growth




Figure 7. Valita® STABILITY workflow. Clones are monitored using ChemStress® before or during the early stages of stability trials to assess whether they are extreme and therefore likely to be unstable.

rate drops but its protein production rate increases, or vice versa, in such a way that their combined effect on titer change is negligible. This 'apparent titer stability', which is masking real underlying biological change, is abundant in our data. Quixotic targets make for challenging prediction problems. An alternative, cleaner formulation of the stability prediction challenge, still aimed at streamlining trials, is to detect clones achieving high titers predominantly through very high growth rates or very high protein production rates as our data indicate such extreme strategies are difficult to maintain over time. By contrast, we observe that high titer clones built from 'moderate' growth and protein production rates have a higher chance of stability than these 'extreme' clones. Our algorithms leverage ChemStress® to detect extreme clones with a high probability of being unstable.




Figure 8. Valita® STABILITY Scoring. Growth and protein production rates for the most extreme clone relative to background distributions for all clones in the same panel, depicted as box-and-whisker plots.





Alongside using ChemStress® to predict industry's established stability metrics, ChemStress® itself can be used to define a rich, functionally relevant measure of clone stability. A clone can be scored for any stress chemical on the ChemStress® plate by calculating how far removed the clone's measurement is from the median measurement derived from a background set of clones. A clone that is many deviations removed from the median is, by definition, extreme relative to the background with respect to that particular stress chemical. By summing scores across all stressors, it is possible to rank clones by relative stability.

Predicting stability metrics is challenging as there is considerable variability in how trials are conducted and how stability is defined, with 'apparent titer stability' particularly clouding the issue. Deploying ChemStress® to inform 'extreme' clone prediction can help streamline trials. Scoring 'extreme' clones directly using ChemStress® provides an additional, function-focused stability metric with the option of dissecting the rich fingerprint to scrutinise the underlying biological mechanisms driving instability.

### Biomanufacturing Process Control

Exploiting the cellular responses to stressful culture environments is a novel and information-rich approach to addressing the key issues of bioprocessing. The inherent variability of cells is a critical component of the biopharmaceutical industry, which can be exploited to create information-rich data sets. New technologies targeted at decreasing the time and costs associated with the development of biotherapeutics must ultimately be prepared to understand the molecular mechanisms of the cells at the centre of bioprocessing. This is the central approach taken by the ChemStress® platform. The assay generates highly relevant information on how a clone reacts to biological stressors, which can be used to discriminate the variation between clones, the deviation of clones from a stable state and variations in the quality of media environments. This is an assay which offers a rapid and easy way to compile a significant amount of functionally relevant information, which is currently lacking in the biotherapeutic workflow. This kind of strategy is necessary to streamline the CLD process and ultimately decrease the substantial cost and time taken to bring a biotherapeutic to market.

### Acknowledgements

The authors would like to acknowledge support from the Horizon2020 Phase II SME Instrument – BIOSIM 778586 grant.

### REFERENCES

- O'Callaghan PM, McLeod J, Pybus LP, Lovelady CS, Wilkinson SJ, Racher AJ, et al. Cell line-specific control of recombinant monoclonal antibody production by CHO cells. *Biotechnol Bioeng*. 2010 Aug 15;106(6):938–51.
- Kim M, O'Callaghan PM, Droms KA, James DC. A mechanistic understanding of production instability in CHO cell lines expressing recombinant monoclonal antibodies. *Biotechnol Bioeng*. 2011 Oct;108(10):2434–46.
- Grainger RK, James DC. CHO cell line specific prediction and control of recombinant monoclonal antibody N-glycosylation. *Biotechnol Bioeng*. 2013 Nov;110(11):2970–83.
- Fernandez-Martell A, Johari YB, James DC. Metabolic phenotyping of CHO cells varying in cellular biomass accumulation and maintenance during fed-batch culture. *Biotechnol Bioeng*. 2018 Mar;115(3):645–60.
- Davies SL, Lovelady CS, Grainger RK, Racher AJ, Young RJ, James DC. Functional heterogeneity and heritability in CHO cell populations. *Biotechnol Bioeng*. 2013 Jan;110(1):260–74.
- Westmann D. Cost and Impact of Bioburden incident. *Ge Healthcare: GE Healthcare*; 2017 Apr.
- Thompson B, Clifford J, Jenns M, Smith A, Field R, Nayyar K, et al. High-throughput quantitation of Fc-containing recombinant proteins in cell culture supernatant by fluorescence polarization spectroscopy. *Anal Biochem*. 2017 01;534:49–55.
- Nicklin P, Bergman P, Zhang B, Triantafellow E, Wang H, Nyfeler B, et al. Bidirectional transport of amino acids regulates mTOR and autophagy. *Cell*. 2009 Feb 6;136(3):521–34.
- Zhong D, Liu X, Schafer-Hales K, Marcus AI, Khuri FR, Sun S-Y, et al. 2-Deoxyglucose induces Akt phosphorylation via a mechanism independent of LKB1/AMP-activated protein kinase signaling activation or glycolysis inhibition. *Molecular Cancer Therapeutics*. 2008 Apr 1;7(4):809–17.
- Roberts DJ, Tan-Sah VP, Ding EY, Smith JM, Miyamoto S. Hexokinase-II Positively Regulates Glucose Starvation-Induced Autophagy through TORC1 Inhibition. *Molecular Cell*. 2014 Feb;53(4):521–33.
- Spitz DR, Phillips JW, Adams DT, Sherman CM, Deen DF, Li GC. Cellular resistance to oxidative stress is accompanied by resistance to cisplatin: The significance of increased catalase activity and total glutathione in hydrogen peroxide-resistant fibroblasts. *Journal of Cellular Physiology*. 1993 Jul;156(1):72–9.
- Poirier F, Bourin P, Bladier D, Joubert-Caron R, Caron M. Effect of 5-azacytidine and galectin-1 on growth and differentiation of the human b lymphoma cell line bl36. *Cancer Cell Int*. 2001 Dec 17;1(1):2.
- Yun CY, Liu S, Lim SF, Wang T, Chung BYF, Jiat Teo J, et al. Specific inhibition of caspase-8 and -9 in CHO cells enhances cell viability in batch and fed-batch cultures. *Metabolic Engineering*. 2007 Sep;9(5–6):406–18.
- Schuerwegh AJ, Stevens WJ, Bridts CH, De Clerck LS. Evaluation of monensin and brefeldin A for flow cytometric determination of interleukin-1 beta, interleukin-6, and tumor necrosis factor-alpha in monocytes. *Cytometry*. 2001 Jun 15;46(3):172–6.
- Vicetti Miguel RD, Maryak SA, Cherpes TL, Brefeldin A, but not monensin, enables flow cytometric detection of interleukin-4 within peripheral T cells responding to ex vivo stimulation with *Chlamydia trachomatis*. *J Immunol Methods*. 2012 Oct 31;384(1–2):191–5.
- Moon JL, Kim SY, Shin SW, Park J-W. Regulation of brefeldin A-induced ER stress and apoptosis by mitochondrial NADP+-dependent isocitrate dehydrogenase. *Biochemical and Biophysical Research Communications*. 2012 Jan;417(2):760–4.
- Kano F, Sako Y, Tagaya M, Yanagida T, Murata M. Reconstitution of brefeldin A-induced golgi tubulation and fusion with the endoplasmic reticulum in semi-intact chinese hamster ovary cells. *Mol Biol Cell*. 2000 Sep;11(9):3073–87.
- Thompson J, Hipwell E, Loo HV, Bannigan J. Effects of cadmium on cell death and cell proliferation in chick embryos. *Reprod Toxicol*. 2005 Dec;20(4):539–48.
- Lag M, Westly S, Lerstad T, Bjørnsrud C, Refsnes M, Schwarze PE. Cadmium-induced apoptosis of primary epithelial lung cells: involvement of Bax and p53, but not of oxidative stress. *Cell Biol Toxicol*. 2002;18(1):29–42.
- Mitachi T, Mezaki M, Yamashita K, Itagaki H. Acidic conditions induce the suppression of CD86 and CD54 expression in THP-1 cells. *J Toxicol Sci*. 2018;43(5):299–309.
- Zou W, Yan M, Xu W, Huo H, Sun L, Zheng Z, et al. Cobalt chloride induces PC12 cells apoptosis through reactive oxygen species and accompanied by AP-1 activation. *J Neurosci Res*. 2001 Jun 15;64(6):646–53.
- Adler-Abramovich L, Vaks L, Carny O, Trudler D, Magno A, Caflisch A, et al. Phenylalanine assembly into toxic fibrils suggests amyloid



- etiology in phenylketonuria. *Nature Chemical Biology*. 2012 Aug;8(8):701–6.
23. Tagde A, Singh H, Kang MH, Reynolds CP. The glutathione synthesis inhibitor buthionine sulfoximine synergistically enhanced melphalan activity against preclinical models of multiple myeloma. *Blood Cancer Journal*. 2014 Jul;4(7):e229–e229.
  24. Koczor CA, Shokolenko IN, Boyd AK, Balk SP, Wilson GL, LeDoux SP. Mitochondrial DNA Damage Initiates a Cell Cycle Arrest by a Chk2-associated Mechanism in Mammalian Cells. *Journal of Biological Chemistry*. 2009 Dec 25;284(52):36191–201.
  25. Hall MD, Marshall TS, Kwit ADT, Miller Jenkins LM, Dulcey AE, Madigan JP, et al. Inhibition of Glutathione Peroxidase Mediates the Collateral Sensitivity of Multidrug-resistant Cells to Tiopronin. *Journal of Biological Chemistry*. 2014 Aug 1;289(31):21473–89.
  26. Li L, Sun Y, Liu J, Wu X, Chen L, Ma L, et al. Histone deacetylase inhibitor sodium butyrate suppresses DNA double strand break repair induced by etoposide more effectively in MCF-7 cells than in HEK293 cells. *BMC Biochemistry*. 2015;16(1):2.
  27. Shen D, Kiehl TR, Khattak SF, Li ZJ, He A, Kayne PS, et al. Transcriptomic responses to sodium chloride-induced osmotic stress: A study of industrial fed-batch CHO cell cultures. *Biotechnology Progress*. 2010;NA-NA.
  28. Lao M-S, Toth D. Effects of Ammonium and Lactate on Growth and Metabolism of a Recombinant Chinese Hamster Ovary Cell Culture. *Biotechnology Progress*. 1997 Oct 7;13(5):688–91.
  29. Zhao Z, Han F, Yang S, Wu J, Zhan W. Oxamate-mediated inhibition of lactate dehydrogenase induces protective autophagy in gastric cancer cells: Involvement of the Akt–mTOR signaling pathway. *Cancer Letters*. 2015 Mar;358(1):17–26.
  30. Morita A, Zhu J, Suzuki N, Enomoto A, Matsumoto Y, Tomita M, et al. Sodium orthovanadate suppresses DNA damage-induced caspase activation and apoptosis by inactivating p53. *Cell Death and Differentiation*. 2006 Mar;13(3):499–511.
  31. Wulhfard S, Baldi L, Hacker DL, Wurm F. Valproic acid enhances recombinant mRNA and protein levels in transiently transfected Chinese hamster ovary cells. *Journal of Biotechnology*. 2010 Jul;148(2–3):128–32.
  32. Yang WC, Lu J, Nguyen NB, Zhang A, Healy NV, Kshirsagar R, et al. Addition of Valproic Acid to CHO Cell Fed-Batch Cultures Improves Monoclonal Antibody Titers. *Molecular Biotechnology*. 2014 May;56(5):421–8.
  33. Curriden S, Englesberg E. Inhibition of growth of proline-requiring Chinese hamster ovary cells (CHO-k1) resulting from antagonism by a system amino acids. *J Cell Physiol*. 1981 Feb;106(2):245–52.
  34. Krampe B, Al-Rubeai M. Cell death in mammalian cell culture: molecular mechanisms and cell line engineering strategies. *Cytotechnology*. 2010 Jun;62(3):175–88.
  35. Codogno P, Meijer AJ. Autophagy and signaling: their role in cell survival and cell death. *Cell Death and Differentiation*. 2005 Nov;12:1509–18.
  36. Calvet A, Li B, Ryder AG. A rapid fluorescence based method for the quantitative analysis of cell culture media photo-degradation. *Analytica Chimica Acta*. 2014 Jan;807:111–9.
  37. Zang L, Frenkel R, Simeone J, Lanan M, Byers M, Lyubarskaya Y. Metabolomics Profiling of Cell Culture Media Leading to the Identification of Riboflavin Photosensitized Degradation of Tryptophan Causing Slow Growth in Cell Culture. *Analytical Chemistry*. 2011 Jul;83(13):5422–30.



## Dr Karen Coss

Dr Karen Coss is a Product Manager at Valitacell. Prior to joining Valitacell, Karen was a Postdoctoral Fellow at King's College London, in the Department of Immunology, Infection and Inflammatory Disease, studying the glycosylation of HIV-1 with an application towards vaccine design. She has also used genomic and glycomic approaches to develop biomarkers to monitor patients with rare metabolic disorders. Karen holds a PhD from University College Dublin, Ireland.

Email: [karen.coss@valitacell.com](mailto:karen.coss@valitacell.com)



## Dr Ben Thompson

Dr Ben Thompson is the Head of Informatics at Valitacell. He is an experienced cell line development scientist with expertise in bioprocess engineering and data science. He has been centrally involved in all Valitacell product innovations since inception and is a co-inventor of Valitacell IP. His interests include fluorescence polarisation technologies for assay development and machine learning for bioprocess optimisation. Ben holds a PhD from the University of Sheffield, UK.

Email: [ben.thompson@valitacell.com](mailto:ben.thompson@valitacell.com)



## Dr Paul Dobson

Dr Paul Dobson is a Data Analyst with Valitacell. Before working with Valitacell he completed numerous postdoctoral positions covering a variety of key subject areas relevant to bioengineering, including bioinformatics, cheminformatics and systems biology. In addition to this, he was a lecturer in bioengineering at the University of Sheffield. Paul holds a PhD in Computational Biology with Machine Learning from the University of Manchester, UK.

Email: [paul.dobson@valitacell.com](mailto:paul.dobson@valitacell.com)



## Dr Jerry Clifford

Dr Jerry Clifford is the co-founder and COO of Valitacell. Prior to Valitacell, Jerry was the commercial director of the Technopath group (IVD quality control manufacturer). He was also Vice President for development at the Institute of Technology Tralee, Ireland. Jerry holds a PhD from Royal College of Surgeons, Ireland and an MBA from University College Cork, Ireland.

Email: [jerry.clifford@valitacell.com](mailto:jerry.clifford@valitacell.com)



## Prof David James

Prof David James is Professor of Bioprocess Engineering in the Department of Chemical and Biological Engineering and the Director of the Advanced Biomanufacturing Centre at the University of Sheffield, UK. His current research is focused on the functional genomics of CHO cell factory function, cellular process design and engineering biopharmaceutical production using computational modelling, systems and synthetic biology. David holds a PhD from King's College London, UK.

Email: [d.c.james@sheffield.ac.uk](mailto:d.c.james@sheffield.ac.uk)